

CASUALTY SIMULATION KIT



Functional Description:

The Casualty Simulation Kit consists of 100 stick-on wound moulages and a variety of makeup material to be applied, in accordance with the accompanying instructions, to selected persons who act as casualties during simulation of a battle or disaster.

After application of the selected wound moulages, makeup is applied to represent different types of bleeding (depending on the type or severity of the wound), frothing at the mouth, loss of stomach content, shock, perspiration, bruises, contusion, and other effects.

The Instructors Guide for the Casualty Simulation Kit provides many full-color illustrations which amplify the text and greatly assist the person applying the moulages and makeup. This Guide lists all makeup and other materials required, along with commercial source of national stock number.

Purpose of Trainer:

The device provides materials for realistically simulating many types of wounds incurred in battle or disaster. The materials are used in training military and civilian groups in first aid to casualties and to condition these groups not to be shocked at the appearance of the casualties. The specific training requirements supported are shown following the descriptive data.

Physical Information:

Kit: 15" x 9" x 14"; 25 lb.

Equipment Required, Not Supplied:

As listed in Applicable Publications.

Special Installation Requirements:

None

Power Requirements:

None

DVC 08-14
NSN 6910-01-300-6541

Training Category/Level Utilized:

Medical/Level 3

Training Requirements Supported:

AMEDD ARTEP with patient care mission

MOSC 91 CMF

SM 081-XXX Tasks

1002 1478 2101 2102
1467

Combat Life Saver Program (All levels)

Applicable Publications:

NAVEXOS P-2709, Instructors Guide for Casualty Simulation Kit, device 11E10

Reference Publications:

FM 8-23

Source and Method of Obtaining:

Available through local TSC

Logistic Responsible Command, Service, or Agency:

USAMMA