

MEDICAL

8
SERIES


- **DEVICE**
- **GTAs**
- **VIDEOTAPES**



<http://pao.hood.army.mil/TSC>

DOLL, NURSING TECHNIQUE TRAINING

DVC 08-01



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
None

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
Intended to be used by military
medical personnel at all levels of
medical training and evaluation of
basic nursing care.

POWER REQUIREMENTS:
None

PHYSICAL INFORMATION:
Life size adult.

PUBLICATIONS:
Manufacturer's
Instruction Guide

WAR-WOUND MOULAGE SET

DVC 08-04



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
None

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
This device is an enhancement of device "Advance Military Make-Up Casualty Simulation Kit" and provides the means to train military personnel in first aid treatment of combat-type casualties and permits trainees to practice and develop first aid skills.

**POWER
REQUIREMENTS:**
None

PHYSICAL INFORMATION:
Life size - 29 lbs - 2.2 cu. ft.

PUBLICATIONS:
AHS TTP 8-3 Pamphlet
"How to Moulage"

INTRAVENOUS THERAPY TRAINER

DVC 08-05
11E1



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
Needle

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
To provide safe, realistic, and
painless practice in performing
intravenous injections in the human
antecubital area.

POWER REQUIREMENTS:
None

PHYSICAL INFORMATION:
Dimensions approximately 24" x 6"

PUBLICATIONS:
None

HYPODERMIC SYRINGE TRAINER

**DVC 08-08
11E4**



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
None

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
The trainer demonstrates the characteristics and operation of the standard 10cc hypodermic syringe, including demonstration of its assembly and of the operation of its component parts.


POWER REQUIREMENTS:
None

PUBLICATIONS:
None

PHYSICAL INFORMATION:
Syringe: Approximately 24" long; 4 lbs.
Container: 55 lbs.

ORAL THERMOMETER TRAINER

**DVC 08-09
11E5**



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
None

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
**For general medical training
demonstrating techniques of
handling and reading the oral
clinical thermometer.**

POWER REQUIREMENTS:
None

PHYSICAL INFORMATION:
30" x 4" x 4"; 12 lbs.

PUBLICATIONS:
None

BUTTOCKS MOULAGE KIT

**DVC 08-11
11E7**



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
None

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
The device permits medical personnel trainees to practice in the intramuscular administration of drugs.

**POWER
REQUIREMENTS:**
None

PHYSICAL INFORMATION:
Life size

PUBLICATIONS:
NAVEXOS P-1946

CASUALTY SIMULATION KIT

DVC 08-14
11E10



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
None

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
The device provides materials for realistically simulating many types of wounds incurred in battle or disaster.

**POWER
REQUIREMENTS:**
None

PHYSICAL INFORMATION:
Kit 15" x 9" x 14"; 25 lbs.

PUBLICATIONS:
FM 8-23

MANIKIN, RESUSCITATION TRAINING

DVC 08-15



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

PURPOSE OF TRAINER:
The device provides realistic practice in mouth-to-mouth and/or cardiopulmonary resuscitation.

PHYSICAL INFORMATION:
Life size, non-inflatable body with removable head.

EQUIPMENT REQUIRED, SUPPLIED:

Letter
"Recommendations for Decontaminating Manikins Used in CPR Available at local TSC Training, 1983 Update."

SPECIAL INSTALLATION REQUIREMENTS:

None

POWER REQUIREMENTS:

10 AA's or 6 D Cell depending on model provided.

PUBLICATIONS:

Commercial Manual

MANIKIN, RESUSCITATION TORSO

DVC 08-16



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

SOURCE AND METHOD OF OBTAINING:
Available at local TSC

PURPOSE OF TRAINER:
To provide a realistic human blood circulation system to be used to demonstrate and train both artificial ventilation and external cardiac compressions.

PHYSICAL INFORMATION:
Device: 24" x 15" x 8" Weight 35 lbs
Carrying Case: 26" x 17" x 10"

EQUIPMENT REQUIRED, SUPPLIED:
Letter
"Recommendations for Decontaminating Manikins Used in CPR Training, 1983 Update."

SPECIAL INSTALLATION REQUIREMENTS:
None

POWER REQUIREMENTS:
None

PUBLICATIONS:
FM 8-230, FM 21-11
STP-21-1-SMCT

MANIKIN, RESUSCITATION TRAINING INFANT

DVC 08-20



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

SOURCE AND METHOD OF OBTAINING:
Available at local TSC

PURPOSE OF TRAINER:
The training device provides realistic practice in mouth-to-mouth and/or cardiopulmonary resuscitation on infants.

PHYSICAL INFORMATION:
Infant size, 22" x 11" x 6",
noninflatable body with carrying case;
weight: 11 lbs.

**EQUIPMENT REQUIRED,
SUPPLIED:**

Letter "Recommendations
for Decontaminating
Manikins Used in CPR
Training, 1983 Update."
Battery, 1.5 V DC (size
"AA") 10 each.

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

POWER REQUIREMENTS:
Battery, 1.5 V DC
(size "AA") 10 each.

PUBLICATIONS:
FM 8-230, FM 21-11,
STP-21-1-SMCT

ARRHYTHMIA SIMULATOR

DVC 08-21



TRAINING CATEGORY/LEVEL UTILIZED:
Medical/Level 3

SOURCE AND METHOD OF OBTAINING:
Available at Local TSC

PURPOSE OF TRAINER:

The device will be used to train field medical care personnel in the recognition of life-threatening cardiac arrhythmias through the use of the Defibrillator/Electrocardiograph/Monitor Recorder NSN: 6515-01-173-2053.

PHYSICAL INFORMATION:

Weight 12 ounces.
Is 1 1/2" x 5 1/2" x 7 1/2" and has 3 universal ECG jacks (for arrhythmias) 10 banana jacks (for ECG testing)

EQUIPMENT REQUIRED, NOT SUPPLIED:

Defibrillator/ECG/
Monitor Recorder,
NSN 6515-01-173-2053

SPECIAL INSTALLATION REQUIREMENTS:

None

POWER REQUIREMENTS:

9 V DC for the device,
110-220 V AC for
defibrillator

PUBLICATIONS:

Operator Manual and
Photo of Abnormal
Tracings

TRAINING KIT, NERVE AGENT ANTIDOTE, MARK I, TYPE I

DVC 08-36



TRAINING CATEGORY/LEVEL UTILIZED:
General/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
None

SOURCE AND METHOD OF OBTAINING:
Available at local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
The nerve agent antidote auto injector training kit is for training individuals in the proper administration of antidote for Chemical Nerve Agent poisoning without using real injectors.

**POWER
REQUIREMENTS:**
None

PHYSICAL INFORMATION:
Each auto-injector incorporates a removable safety cap. The assembled kit comes packaged in a Styrofoam case.

PUBLICATIONS:
FM 21-2-SMCT
STP 21-1-SMCT

TRAINING DEVICE, DIAZEPAM AUTO INJECTION

DVC 08-37



TRAINING CATEGORY/LEVEL UTILIZED:
General/Level 3

**EQUIPMENT REQUIRED,
NOT SUPPLIED:**
None

SOURCE AND METHOD OF OBTAINING:
Available at local TSC

**SPECIAL INSTALLATION
REQUIREMENTS:**
None

PURPOSE OF TRAINER:
The nerve agent antidote auto injector training kit is for training individuals in the proper administration of antidote for Chemical Nerve Agent poisoning without using real injectors. The Diazepam Auto injector is for buddy-aid.

**POWER
REQUIREMENTS:**
None

PUBLICATIONS:
FM 21-2-SMCT
STP 21-1-SMCT

PHYSICAL INFORMATION:
Each device incorporates a removable safety cap. Kit is contained in a plastic case.